

POWER MOBILITY WHEELCHAIR

ESCAPE HP1

A steel-framed model serves as an economical alternative.



ESCAPE DX HP5

An Aluminum-framed model features Aluminum die-casting lugs instead of welded ones.



ESCAPE LX HP8



Adjustable Backrest
(90~110 degrees)

HP8 Escape LX is a strong, durable electric wheelchair ideal for portability and general indoor/outdoor use. Available in 16", 18", and 20" seat widths. Adjustable backrest



ESCAPE S P12S

A model incorporates a rear-wheel suspension system, and thus provides great comfort and ensures a good ride.



ESCAPE SX P12SX

A strong and enduring model, with 9" front castors and 14" rear wheels installed, is ready for a faster stroll.



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
HP1	110kg	1090mm x 630mm	62kg	DC24V / 200W	6KPH	12V 36Ah x 2pcs	30km	12.5" x 2.5"	5Amp Off Board	820mm
HP5	135kg	1100mm x 685mm	62kg	DC24V / 200W	6KPH	12V 36Ah/50Ah x 2pcs	30km	12.5" x 2.5"	5Amp Off Board	820mm
HP8	135kg	1085mm x 640mm	66kg	DC24V / 200W	6KPH	12V 36Ah/50Ah x 2pcs	30km/45km	12.5" x 2.5"	5Amp Off Board	820mm

Heartway reserves the right to modify or change the specifications without notice.

Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P12S	135kg	1120mm x 610mm	73kg	DC 24V / 200W	6KPH	12V 60Ah x 2pcs	45km	12.5" x 2.5"	5Amp Off Board	1240mm
P12SX	160kg	1130mm x 660mm	83kg	DC 24V / 600W	12KPH	12V 50Ah x 2pcs	45km	14" x 3"	5Amp Off Board	1240mm

Heartway reserves the right to modify or change the specifications without notice.

POWER MOBILITY WHEELCHAIR



ESCAPE HD HP1HD



A heavy-duty steel-framed wheelchair, structured with a double-cross brace design and equipped with 4-pole high power motors, offers seat width options of 22" and 24" for the heavyweights.

ESCAPE LXHD HP8HD



HP8HD is a new addition to Heartway's renown Bariatric range of products. Its double-cross bar design greatly strengthens the structure of the product. It offers seat width options of 24", 26" and 28". Also as a folding chair, at the same time it offers high portability for our big friends who don't mind travelling a little bit!



Specification	Weight Capacity	Dimension LxW	Weight W/Battery	Motor Type	Max speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
HP1HD	180kg	1150mm x 750mm	82kg	DC 24V / 700W	10KPH	12V 50Ah x 2pcs	30km	14" x 3"	5Amp Off Board	1000mm
HP8HD	225kg	1170mm x 910mm	86kg	DC 24V / 700W	10KPH	12V 50Ah x 2pcs	30km	14" x 3"	8Amp Off Board	1140mm

MAXX P3D



A mid-wheel drive power wheelchair, coming with a unique independent suspension system, offers great stability and lots of comfort and is featured with an obstacle height of 6cm and more.



MAXX R P3DR



A mid-wheel drive power wheelchair with a rehab seating system provides a wonderful driving experience due to its very special independent suspension system.



MAXX RT P3DRT



A mid-wheel drive model featuring a tilt-in-space rehab seating system offers great mobility not just with the chair but also on the chair.



Specification	Weight Capacity	Dimension LxW	Weight W/Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P3D	160kg	1080mm x 680mm	107kg	DC 24V / 450W	8KPH	12V 50Ah x 2pcs	30km	12" x 4"	5Amp Off Board	550mm
P3DR	160kg	1080mm x 680mm	118kg	DC 24V / 450W	8KPH	12V 50Ah x 2pcs	30km	12" x 4"	5Amp Off Board	550mm
P3DRT	160kg	1080mm x 680mm	133kg	DC 24V / 450W	8KPH	12V 50Ah x 2pcs	30km	12" x 4"	5Amp Off Board	550mm

Heartway reserves the right to modify or change the specifications without notice.