

POWER MOBILITY WHEELCHAIR

P35

- Lightweight travel chair weighs only 29kg .
- One step easy folding- perfect for loading/unloading in vehicle.
- Anti-theft battery covers to protect the battery from being stolen.
- Takes up little space for storage once folded.
- Driving range of 20 km. Weight capacity of 115 kg



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
P35	115kg	780mm x 590mm	29kg	DC24 / 180W	8KPH	8 AH (Lithium)	20km	9" x 2.5"	2Amp Off Board	1390mm

Heartway reserves the right to modify or change the specifications without notice.