

# POWER MOBILITY SCOOTER



## PF7S

### ROYALE 4S



Fully independent suspensions and a top speed up to 16KPH  
Equipped with a four-pole 1500W motor and a 200 Amp controller  
Heavy duty battery upgrade available to increase travelling range  
Digital display with odometer, speed, trip distance, battery indicator  
Double seat available to accommodate more than one passenger  
Equipped with G-Sensor for speed reduction at corner  
High quality sunproof canopy with luxury design



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
PF7S	160kg	1690mm x 750mm	211kg	DC24V 1500W	16KPH	12V 100Ah x 2 pcs	60km	19" x 4.7"	8 Amp Off Board	2500mm