

BRIO S19



Folds in a simple movement

Lightweight design. Easy to lift into the car boot


Takes up very little space for storage once folded

Easy to lug around at home or shopping malls

Designed to install and uninstall the lithium batteries quickly

Anti-theft battery covers to protect the battery from being stolen

Carries up to 115kg for 15-20 km travelling distance per charge



BRIO S19⁺



The Heartway S19⁺ incorporates compact, foldable mobility scooter design along with additional features including larger pneumatic tires, a more powerful controller, free LED light, higher ground clearance and enhanced batteries (14.5 AH) for longer travel range.

It also provides a rider weight capacity accommodating users weighing up to 115kg.



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
S19	115kg	950mm x 500mm	26kg	DC 24V / 270 W	6KPH	11.5Ah(Lithium)	15km	8" x 2"	2.5Amp Off Board	820mm
S19 ⁺	115kg	950mm x 500mm	26.5kg	DC 24V / 270 W	6KPH	14.5Ah(Lithium)	20km	9" x 2.8"	2.5Amp Off Board	820mm

Heartway reserves the right to modify or change the specifications without notice.