

POWER MOBILITY SCOOTER

BRIO

S19V

A truly advanced, compact and transportable design that folds in one simple movement with the push of a foot pedal system. No parts need to be removed for storage and transportation thanks to lightweight lithium battery technology and the easy one-handed effort to collapse the durable aluminium frame. With the push of a foot pedal system, the Heartway S19V folds up in seconds making it simple to store or transport. Innovative safety features such as a sensor to control cornering speeds ensure the Heartway S19V stands out from the crowd

- 1.RFID swipe card to turn on/turn off scooter power
- 2.Folds in a simple movement with the push of a foot pedal
- 3.Total weight less than 27kg (batteries included)
- 4.Easy to lug around at home, airport, cruises, or shopping malls
- 5.Anti-theft battery covers to protect the battery from being stolen
- 6.Compact turning radius
- 7.Simple to store and transport
- 8.FDA, CE, TUV, and ISO standards approval



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
S19V	115KG	970mm x 530mm	27KG	DC 24V / 270W	6KPH	11.5 Ah(Lithium)	15km	9" x 2.5"	2Amp Off Board	820mm

Heartway reserves the rights to amend or change the specifications if required and necessary. There may not be further notification regarding the amendments or changes.