

# POWER MOBILITY SCOOTER

EASY MOVE

S21



- Automatic folding & unfolding with one touch on a remote control
- Powered by 8.7Ah lithium battery
- Lightweight: 19 kg
- Maximum travel range – up to 11 km
- Features a large 115kg weight carrying capacity
- LED headlight; height adjustable tiller; dual hand two way throttle control
- Comes with carrying bag; easy to carry on the travel
- Ideal for cruise ships, and shopping centers or folding away into the boot of a car
- European CE Certified (EN 12184)



Specification	Weight Capacity	Dimension LxW	Weight W/ Battery	Motor Type	Max Speed	Battery	Battery Range	Drive Wheel	Charger Type	Turning Radius
S21	115kg	860mm x 420mm	19kg	DC24 / 180W	6KPH	8.7 AH (Lithium)	11km	8" x 2"	2Amp Off Board	780mm

Heartway reserves the right to modify or change the specifications without notice.